

Using the beneficial effects of Spark Waves combined with the latest bioengineering technology to treat urogenital disorders

urogold100



Treatment with the urogold100 provides a long-term, effective and highly tolerable solution for patients with urogenital diseases

High performance through a high quality product

The urogold100 device is characterized by its compact size and modern design. This leads to an especially easy transport and comfortable handling of the device during treatment.

Premium materials and the latest MTS technology guarantee highquality Spark Waves for an ideal treatment process.

MTS offers different applicators for the various indications ranging from an unfocused applicator (OP155) for the treatment of CPPS and ED to an either soft-focused (OE035) or focused applicator (OE050) for the treatment of IPP.

Beneficial features of urogold100

The advantages of the treatment with the urogold100 device are its versatile and professional applicators which make it easy to repeat the treatment on a regular basis. The therapy of urogenital indications is uncomplicated to conduct and highly tolerable for the patient. Spark Wave therapy features a long-term effectiveness so as to provide a high degree of satisfaction for both patients and physicians.





Making use of an established technology

Spark Wave Therapy has been proven in orthopaedics and wound treatment for many years now. Through extensive research, the beneficial effects of Spark Waves have gradually been extended to the urogenital field. With their latest innovation urogold100 MTS offers a gentle and effective long-term treatment method for patients with urogenital disorders.

Indications for urogold100

- Erectile dysfunction (ED)
- Peyronie's disease IPP (Induratio penis plastica)
- Chronic pelvic pain syndrome (CPPS) / Chronic prostatitis

The bioengineering aspect of urogold100

The urogold100 represents an innovative treatment option that utilizes the biological effects of Spark Waves by stimulating the self-healing mechanism of the affected tissue on cellular level. Spark Waves promote the generation of new blood vessels and the release of angiogenic growth factors. It has been shown that they also initiate the migration of stem cells and therefore help improve the blood flow in the treated area. In addition, Spark Wave-induced revascularization processes can alleviate pain. The combination of all mentioned factors leads to reductions in muscle tension, spasticity and pain for CPPS as well as the long-term recovery of the erection function of the patient. The reduction of plague size and improvement of penile angulation is another healing factor for patients with IPP.

Comfortable treatment on an outpatient basis

The treatment of IPP, ED and CPPS is carried out on an outpatient basis. Dependent on the indication, the Spark Waves are applied onto the affected area that has been determined upfront through ultrasound, palpation or anatomical knowledge. Most importantly, the Spark Wave therapy with urogold100 is pain-free and thus requires no anaesthesia. No side effects have been reported for the Spark Wave treatment of urogenital disorders.

Issue March 2013





Technical Data

Applicator OE035

0 - 79 mm Therapeutic Focus fz (-6dB): 0.01 - 0.16 mJ/mm² Energy Flux Density (ED): Focus Energy (E total / -6dB): 0.34 - 4.22 mJ

Applicator OE050

Therapeutic Focus fz (-6dB): 0 - 62 mm 0.03 - 0.27 mJ/mm² Energy Flux Density (ED): 0.30 - 6.15 mJ Focus Energy (E total / -6dB):

Applicator OP155

Therapeutic Focus fz (-6dB): 0 - 82 mm 0.01 - 0.19 mJ/mm² Energy Flux Density (ED): Focus Energy (E total / -6dB): 0.40 - 3.88 mJ

Dimensions: 21.8 x 40 x 45.9 cm Power Supply: AC 100 - 240 V, 50 - 60 Hz 200 VA Power Input: 0.5 - 8 Hz Repetition Rate: Weight: 16.5 kg



